

Antibodies Product Datasheet

Product Name:	Keratin(50-60kDa) mouse mAb
Catalog Number:	N0974
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb dilution 1:2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human KRT5 beta protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	krt5
Protein Name:	
Other Name:	58 kDa cytokeratin; CK 5; CK-5; CK5; Cytokeratin-5; Cytokeratin5; DDD; DDD1; EBS 2; EBS2; epidermolysis bullosa simplex 2 Dowling-Meara/Kobner/Weber-Cockayne types; K2C5_HUMAN; K5; Keratin 5 (epidermolysis bullosa simplex Dowling-Meara/Kobner/Weber-Cockayne)
Human Swiss Prot No.:	P13647
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Keratin(50-60kDa) proteins.
Storage and Stability:	-20°C/1 year