

Antibodies Product Datasheet

Product Name:	Lamin A mouse mAb
Catalog Number:	N0975
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human LMNA protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	Imna
Protein Name:	
Other Name:	70 kDa lamin ; CDCD1 ; CDDC ; CMD1A ; CMT2B1 ; EMD2 ; FPL ; FPLD ; HGPS ; IDC ; LAMIN A ; lamin A/C ; LAMIN C ; Lamin-A/C ; LDP1 ; LFP ; LGMD1B ; LMN 1 ; LMN A ; LMN C ; LMN1 ; LMNA ; LMNA_HUMAN ; LMNC ; NY REN 32 antigen ; PRO1 ; Renal carcinoma antigen
Human Swiss Prot No.:	P02545
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Lamin A and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year