

Antibodies Product Datasheet

Product Name:	BRG1 mouse mAb
Catalog Number:	N0977
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human BRG1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	smarca4
Protein Name:	
Other Name:	ATP dependent helicase SMARCA4; ATP-dependent helicase SMARCA4; BAF 190; BAF190; BAF190A; Brahma protein homolog 1; Brahma protein like 1; BRG1; BRG1 associated factor 190A; BRG1 protein; BRG1-associated factor 190A; BRM/SWI2 related gene 1; Global transcription acti
Human Swiss Prot No.:	P51532
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRG1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year