

Antibodies Product Datasheet

Product Name:	FEN-1 mouse mAb
Catalog Number:	N0979
Reactivity:	Human, Mouse, Rat
Applications:	WB, IP
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human FEN-1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	fen1
Protein Name:	
Other Name:	DNase IV; FEN-1; FEN1; FEN1_HUMAN; Flap endonuclease 1; Flap structure specific endonuclease 1; Flap structure-specific endonuclease 1; hFEN-1; hFEN1; Maturation factor 1; MF1; Rad2.
Human Swiss Prot No.:	P39748
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FEN-1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year