

Antibodies Product Datasheet

Product Name:	XRCC1 mouse mAb
Catalog Number:	N0982
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human XRCC1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	xrcc1
Protein Name:	
Other Name:	DNA repair protein XRCC1;RCC;X ray repair complementing defective repair in chinese hamster;X ray repair complementing defective repair in chinese hamster cells 1;X ray Repair Complementing Defective Repair in Chinese Hamster Cells;X ray repair cross comp
Human Swiss Prot No.:	P18887
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of XRCC1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year