

Antibodies Product Datasheet

Product Name:	Hsp60 mouse mAb
Catalog Number:	N0984
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:100 ip dilution 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hsp60 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hsp60
Protein Name:	
Other Name:	60 kDa chaperonin; 60 kDa heat shock protein mitochondrial; 60 kDa heat shock protein, mitochondrial; CH60_HUMAN; Chaperonin 60; Chaperonin; Chaperonin, 60-KD; CPN 60; CPN60; fa04a05; GROEL; GroEL Homolog; Heat shock 60kD protein 1 chaperonin; Heat shock 60kDa protein
Human Swiss Prot No.:	P10809
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp60 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year