

Antibodies Product Datasheet

Product Name:	GADD45 α mouse mAb
Catalog Number:	N0985
Reactivity:	Transfected
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human GADD45 α protein expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	gadd45a
Protein Name:	
Other Name:	DDIT 1; DDIT-1; DDIT1; DNA damage inducible transcript 1; DNA damage-inducible transcript 1 protein; GA45A_HUMAN; GADD45; GADD45A; Growth arrest and DNA damage inducible 45 alpha; Growth arrest and DNA damage inducible alpha; Growth arrest and DNA damage-
Human Swiss Prot No.:	P24522
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected Only.
Storage and Stability:	-20°C/1 year