

Antibodies Product Datasheet

Product Name:	Bag1 mouse mAb
Catalog Number:	N0986
Reactivity:	
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Bag1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	bag1
Protein Name:	
Other Name:	Bag 1;BAG family molecular chaperone regulator 1;BAG-1;BAG1;BAG1_HUMAN;BCL 2 Associated Athanogene 1;BCL 2 associated athanogene;BCL 2 binding athanogene 1;Bcl-2-associated athanogene 1;BCL2 associated athanogene 1;BCL2 associated athanogene;Glucocorticoid
Human Swiss Prot No.:	Q99933
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected
Storage and Stability:	-20°C/1 year