

Antibodies Product Datasheet

Product Name:	Human IgG mouse mAb
Catalog Number:	N0988
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Human IgG protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	human igg
Protein Name:	
Other Name:	FLJ39988; FLJ40036; FLJ40253; FLJ40587; FLJ40789; FLJ40834; G1m marker; G2m marker; G3m marker; G4m marker; HDC; Heavy chain disease protein; Human immunoglobulin G; Ig gamma 1 chain C region; Ig gamma 2 chain C region; Ig gamma 3 chain C region.
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects the heavy chain of Human IgG and does not cross-react with Human IgM.
Storage and Stability:	-20°C/1 year