

Antibodies Product Datasheet

Product Name:	Rabbit IgG mouse mAb
Catalog Number:	N0989
Reactivity:	Rabbit
Applications:	WB,ELISA
Dilutions:	wb dilution 1:50000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Rabbit IgG
Concentration:	1 mg/ml
Gene Name:	rabbit igg
Protein Name:	
Other Name:	Antibody; Ig gamma chain C region; Ig kappa-b4 chain C region; Ig kappa-b5 chain C region; Ig kappa-b9 chain C region; Ig lambda chain C region; K-BAS.
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects rabbit IgG.
Storage and Stability:	-20°C/1 year