

Antibodies Product Datasheet

Product Name:	H-FABP mouse mAb
Catalog Number:	N1007
Reactivity:	
Applications:	sELISA, Detector
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant H-FABP protein expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	fabp3
Protein Name:	
Other Name:	422 protein; Cardiac Fatty Acid Binding Protein; FABP 11; FABP 3; FABP11; FABP3; FABPH_HUMAN; fatty acid binding protein 11; Fatty acid binding protein 3; Fatty acid binding protein 3 muscle and heart; Fatty acid binding protein 3 muscle and heart mammary
Human Swiss Prot No.:	P05413
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	-20°C/1 year