

Antibodies Product Datasheet

Product Name:	eIF2 α mouse mAb
Catalog Number:	N1009
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, IP
Dilutions:	wb dilution 1:1000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human eIF2 α protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	eif2s1
Protein Name:	
Other Name:	Dmel\CG6376 ; Dmel_CG6376 ; drosE2F1 ; E(Sev-CycE)3A ; E(var)3-93E ; E2-promoter binding facto ; E2F 1 ; E2F transcription factor 1 ; E2F-1 ; E2f-PA ; E2f-PB ; E2f-PC ; E2F1 ; E2f1 E2F transcription factor 1 ; E2F1_HUMAN ; Evar(3)164 ; KIAA4009 ; l(3)07172 ; l(3)j3B1 ; l(3)
Human Swiss Prot No.:	P05198
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of eIF2 α and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year