

Antibodies Product Datasheet

Product Name:	Leptin mouse mAb
Catalog Number:	N1011
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant full length of human Leptin protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	lep
Protein Name:	
Other Name:	FLJ94114; LEP; LEP_HUMAN; LEPD; Leptin (murine obesity homolog); Leptin (obesity homolog, mouse); Leptin; Leptin Murine Obesity Homolog; Leptin Precursor Obesity Factor; OB; Obese Protein; Obese, mouse, homolog of; Obesity; Obesity factor; Obesity factor;
Human Swiss Prot No.:	P41159
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects human Leptin proteins.
Storage and Stability:	-20°C/1 year