

Antibodies Product Datasheet

Product Name:	Human EGF mouse mAb
Catalog Number:	N1015
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant Human EGF protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	egf
Protein Name:	
Other Name:	AI790464 ;Beta urogastrone ;Epidermal Growth Factor ;HOMG4 ;OTTHUMP00000219721 ;OTTHUMP00000219722 ;Pro epidermal growth factor ;Pro epidermal growth factor precursor (EGF) ;URG ;Urogastrone .
Human Swiss Prot No.:	P01133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects Human EGF.
Storage and Stability:	-20°C/1 year