

Antibodies Product Datasheet

Product Name:	4E-BP2 mouse mAb
Catalog Number:	N1019
Reactivity:	Transfected
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein of human 4E-BP2.
Concentration:	1 mg/ml
Gene Name:	eif4ebp2
Protein Name:	
Other Name:	4E-BP2 ; 4EBP2 ; 4EBP2_HUMAN ; eIF4E binding protein 2 ; eIF4E-binding protein 2 ; Eif4ebp2 ; Eukaryotic translation initiation factor 4E binding protein 1 ; Eukaryotic translation initiation factor 4E-binding protein 2 ; PHASII ; phosphorylated, heat and
Human Swiss Prot No.:	Q13542
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected Only.
Storage and Stability:	-20°C/1 year