

Antibodies Product Datasheet

Product Name:	EIF4EBP3 mouse mAb
Catalog Number:	N1020
Reactivity:	Transfected
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein of human EIF4EBP3.
Concentration:	1 mg/ml
Gene Name:	eif4ebp3
Protein Name:	
Other Name:	4E BP3; 4E-BP3; 4EBP3; eIF4E binding protein 3; Eukaryotic initiation factor 4E binding protein 3; Eukaryotic translation initiation factor 4E binding protein 3; OTTHUMP00000159879.
Human Swiss Prot No.:	O60516
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected Only.
Storage and Stability:	-20°C/1 year