

Antibodies Product Datasheet

Product Name:	CK-MM mouse mAb
Catalog Number:	N1022
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CK-MM protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ckmm
Protein Name:	
Other Name:	CKM; CKMM; Creatine kinase M; Creatine kinase M chain; Creatine kinase M type; Creatine kinase M-type; Creatine kinase muscle; Creatine kinase, muscle type; KCRM_HUMAN; M-CK; MCK; Muscle creatine kinase.
Human Swiss Prot No.:	P06732
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects Human CK-MM.
Storage and Stability:	-20°C/1 year