

Antibodies Product Datasheet

Product Name:	MKK3/6 mouse mAb
Catalog Number:	N1025
Reactivity:	Human,Rat
Applications:	WB,IP
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human MKK6 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	map2k3
Protein Name:	
Other Name:	AW212142;dual specificity mitogen activated protein kinase kinase 3;Dual specificity mitogen activated protein kinase kinase 6;Dual specificity mitogen-activated protein kinase kinase 3;MAP kinase kinase 3;MAP kinase kinase 6;map2k3;MAP2K6;MAPK ERK kinase
Human Swiss Prot No.:	P52564
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MKK3 and MKK6 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year