

Antibodies Product Datasheet

Product Name:	Anthrax Protective Antigen mouse mAb
Catalog Number:	N1026
Reactivity:	Bacillus anthracis
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	
Concentration:	1 mg/ml
Gene Name:	pad4
Protein Name:	
Other Name:	anthrax PA ;Anthrax toxins translocating protein ;PA 83 ;PA ;PA83 ;Pag ;PagA ;Protective antigen.
Human Swiss Prot No.:	P13423
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects Anthrax Protective Antigen.
Storage and Stability:	-20°C/1 year