

Antibodies Product Datasheet

Product Name:	Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) mouse mAb
Catalog Number:	N1029
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb dilution 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) synthetic peptide conjugated to KLH.
Concentration:	1 mg/ml
Gene Name:	mapk
Protein Name:	
Other Name:	ERK 1;ERK 2;ERK-2;ERK1;erk1/2;ERK2;ERT1;ERT2;Extracellular signal regulated kinase 1;Extracellular signal regulated kinase 1;Extracellular signal regulated kinase 2;Extracellular signal regulated kinase 2;Extracellular signal-regulated kinase 2;HS44KDAP;H
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204).
Storage and Stability:	-20°C/1 year