

Antibodies Product Datasheet

Product Name:	Hsp90 α mouse mAb
Catalog Number:	N1031
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hsp90 alpha protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	hsp90aa1
Protein Name:	
Other Name:	Heat shock 86 kDa; Heat shock 90kDa protein 1 alpha; Heat shock protein 90kDa alpha (cytosolic) class A member 1; heat shock protein 90kDa alpha (cytosolic), class A member 2; Heat shock protein HSP 90-alpha; HS90A_HUMAN; HSP 86; HSP 86; HSP 86; HSP86; Hs
Human Swiss Prot No.:	P07900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp90 alpha and does not cross-react with Hsp90 beta and other proteins.
Storage and Stability:	-20°C/1 year