

Antibodies Product Datasheet

Product Name:	Keratin8/18 mouse mAb
Catalog Number:	N1033
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	WB 1:500-2000, ICC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein corresponding to full length Keratin8 of human.
Concentration:	1 mg/ml
Gene Name:	krt8/18
Protein Name:	
Other Name:	card2; Cardiac autoantigen 2 120kD; CK 8; CK-8; ck8; Cyk8; cyk8; CYKER; Cytokeratin 8; Cytokeratin endo A; Cytokeratin-8; Cytokeratin8; DreK8; EndoA; k0; k2c8; K2C8_HUMAN; k8; Keratin 8; Keratin type ii cytoskeletal 8; Keratin, type II cytoskeletal 8; Keratin-8; Keratin8;
Human Swiss Prot No.:	P05787/P05783
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of keratin8/18.
Storage and Stability:	-20°C/1 year