

Antibodies Product Datasheet

Product Name:	GOT2 mouse mAb
Catalog Number:	N1035
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein corresponding to fragment (aa/225-430) GOT2 .
Concentration:	1 mg/ml
Gene Name:	got2
Protein Name:	
Other Name:	AATM_HUMAN ; AL022787 ; Aspartate aminotransferase ; Aspartate aminotransferase 2 ; Aspartate aminotransferase, mitochondrial ; ASPATA ; EC 2.6.1.1 ; FABP 1 ; FABP pm ; FABP-1 ; FABPpm ; Fatty acid binding protein ; Fatty acid-binding protein ; FLJ40994 ; Glutamate oxa
Human Swiss Prot No.:	P00505
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GOT2.
Storage and Stability:	-20°C/1 year