

Antibodies Product Datasheet

Product Name:	L-FABP mouse mAb
Catalog Number:	N1040
Reactivity:	Human
Applications:	ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human L-FABP protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	fabp1
Protein Name:	
Other Name:	FABP 1;FABP1;FABP1;FABPL;FABPL_HUMAN;Fatty Acid Binding Protein 1;Fatty acid binding protein 1 liver;Fatty Acid Binding Protein;Fatty acid-binding protein 1;Fatty acid-binding protein;Fatty acid-binding protein liver;L FABP;L-FABP;liver;Liver-type fatty a
Human Swiss Prot No.:	P07148
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	
Storage and Stability:	-20°C/1 year