

Antibodies Product Datasheet

Product Name:	AIF mouse mAb
Catalog Number:	N1044
Reactivity:	Human
Applications:	WB,IF/ICC
Dilutions:	wb dilution 1:1000 icc dilution 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human AIF protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	aifm1
Protein Name:	
Other Name:	AIFM1; AIFM1_HUMAN; Apoptosis inducing factor 1, mitochondrial; Apoptosis inducing factor; Apoptosis inducing factor, mitochondrion associated, 1; Apoptosis-inducing factor 1; CMTX4; COWCK; COXPD6; Harlequin; Hq; mAIF; MGC111425; MGC5706; mitochondrial; N
Human Swiss Prot No.:	O95831
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AIF and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year