

Antibodies Product Datasheet

Product Name:	FADD mouse mAb
Catalog Number:	N1048
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human FADD protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	fadd
Protein Name:	
Other Name:	FADD;FADD protein;FADD_HUMAN;Fas (TNFRSF6) associated via death domain;Fas associated via death domain;Fas associating death domain containing protein;Fas associating death domain containing protein;Fas associating protein;Fas associating protein with dea
Human Swiss Prot No.:	Q13158
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FADD and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year