

Antibodies Product Datasheet

Product Name:	MAP2(N-term) mouse mAb
Catalog Number:	N1052
Reactivity:	Human, Mouse
Applications:	WB, IP, ICC
Dilutions:	wb dilution 1:1000 icc dilution 1:150
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human MAP2 (N-terminus) protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	map2
Protein Name:	
Other Name:	DKFZp686I2148; Dendrite specific MAP; DKFZp686I2148; MAP 2; MAP dendrite specific; MAP-2; MAP2; MAP2A; MAP2B; MAP2C; Microtubule associated protein 2; Microtubule-associated protein 2; Mtap 2; MTAP2_HUMAN.
Human Swiss Prot No.:	P11137
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MAP2 (N-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year