

Antibodies Product Datasheet

Product Name:	Lamin B1 mouse mAb
Catalog Number:	N1055
Reactivity:	Human, Mouse, Monkey
Applications:	WB
Dilutions:	wb dilution 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human Lamin B1 protein.
Concentration:	1 mg/ml
Gene Name:	Lamin B1
Protein Name:	
Other Name:	ADLD; lamin B1; Lamin-B1; LMN; LMN2; LMNB; Lmnb1; LMNB1_HUMAN; MGC111419; OTTHUMP00000159218.
Human Swiss Prot No.:	P20700
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Lamin B1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year