

Antibodies Product Datasheet

Product Name:	NFKB1 mouse mAb
Catalog Number:	N1057
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human NF- κ B1 p105/p50 protein.
Concentration:	1 mg/ml
Gene Name:	NF- κ B1 p105/p50
Protein Name:	
Other Name:	DKFZp686C01211; DNA binding factor KBF1; DNA binding factor KBF1 EBP1; DNA-binding factor KBF1; EBP 1; EBP-1; EBP1; KBF1; MGC54151; NF kappa B; NF kappaB; NF kappabeta; NF kB1; NFkappaB; NFKB 1; NFKB p105; NFKB p50; Nfkb1; NFKB1_HUMAN; Nuclear factor kapp
Human Swiss Prot No.:	P19838
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NF- κ B1 p105/p50 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year