

Antibodies Product Datasheet

Product Name:	Lamin B2 mouse mAb
Catalog Number:	N1058
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb dilution 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human Lamin B2 protein.
Concentration:	1 mg/ml
Gene Name:	Lamin B2
Protein Name:	
Other Name:	LAMB 2; LAMB2; Lamin-B2; LMN 2; LMN B2; LMN2; LMNB 2; LMNB2; LMNB2_HUMAN; MGC2721; RGD1563803.
Human Swiss Prot No.:	Q03252
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Lamin B2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year