

Antibodies Product Datasheet

Product Name:	E-Cadherin mouse mAb
Catalog Number:	N1059
Reactivity:	Human
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	
Concentration:	1 mg/ml
Gene Name:	cdh1
Protein Name:	
Other Name:	Arc 1; CADH1_HUMAN; Cadherin 1; cadherin 1 type 1 E-cadherin; Cadherin1; CAM 120/80; CD 324; CD324; CD324 antigen; cdh1; CDHE; E-Cad/CTF3; E-cadherin; ECAD; Epithelial cadherin; epithelial calcium dependant adhesion protein; LCAM; Liver cell adhesion mole
Human Swiss Prot No.:	P12830
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of E-Cadherin and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year