

Antibodies Product Datasheet

Product Name:	HtrA2/Omi mouse mAb
Catalog Number:	N1061
Reactivity:	Mouse,Hamster
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human HtrA2/Omi protein.
Concentration:	1 mg/ml
Gene Name:	htra2
Protein Name:	
Other Name:	High temperature requirement protein A2;HTRA 2;HtrA like serine protease;HtrA serine peptidase 2;HtrA, E. coli, homolog of, 2;HtrA2;HTRA2_HUMAN;mitochondrial;Omi stress regulated endoprotease;Omi stress-regulated endoprotease;PARK 13;PARK13;Protease serin
Human Swiss Prot No.:	O43464
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HtrA2/Omi and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year