

Antibodies Product Datasheet

Product Name:	COX IV Monoclonal Antibody(6C8), AbFluor™ 488 Conjugated
Catalog Number:	N1067
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICC
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC 1:50-300, IF 1:200 .
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG1
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	COX4I1
Protein Name:	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Other Name:	COX4I1
Human Swiss Prot No.:	P13073
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	COX IV Monoclonal Antibody(6C8) AbFluor™ 488 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.