

Antibodies Product Datasheet

Product Name:	Cyclophilin B Monoclonal Antibody(2B10), AbFluor™ 680 Conjugated
Catalog Number:	N1083
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC, WB, IHC-p
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC 1:50-300, IF: 1:50-200.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	PPIB
Protein Name:	Peptidyl-prolyl cis-trans isomerase B (PPIase B) (EC 5.2.1.8) (CYP-S1) (Cyclophilin B) (Rotamase B) (S-cyclophilin) (SCYLP)
Other Name:	Cyclophilin B
Human Swiss Prot No.:	P23284
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	Cyclophilin B Monoclonal Antibody(2B10) AbFluor™ 680 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.