

Antibodies Product Datasheet

Product Name:	Flag-Tag Monoclonal Antibody(2C5), AbFluor™ 594 Conjugated
Catalog Number:	N1093
Reactivity:	Species independent
Applications:	WB,IP,IF/ICC
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:2000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	Flag Tag
Protein Name:	Flag-Tag DDDDK TAG, DDDDK-TAG
Other Name:	Flag
Human Swiss Prot No.:	
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	Flag-Tag Monoclonal Antibody(2C5) AbFluor™ 594 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.