

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Flag-Tag Monoclonal Antibody(2C5), HRP Conjugated                                                                                                                                                              |
| <b>Catalog Number:</b>        | N1099                                                                                                                                                                                                          |
| <b>Reactivity:</b>            | Species independent                                                                                                                                                                                            |
| <b>Applications:</b>          | WB,IP,IF/ICC                                                                                                                                                                                                   |
| <b>Dilutions:</b>             | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000.                                                                     |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                                                                              |
| <b>Source:</b>                | Monoclonal Mouse IgG                                                                                                                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                                                                                     |
| <b>Immunogen:</b>             |                                                                                                                                                                                                                |
| <b>Concentration:</b>         | 1mg/ml                                                                                                                                                                                                         |
| <b>Gene Name:</b>             | Flag Tag                                                                                                                                                                                                       |
| <b>Protein Name:</b>          | Flag-Tag DDDDK TAG, DDDDK-TAG                                                                                                                                                                                  |
| <b>Other Name:</b>            | Flag                                                                                                                                                                                                           |
| <b>Human Swiss Prot No.:</b>  |                                                                                                                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                                                                                         |
| <b>Specificity:</b>           | Flag-Tag Monoclonal Antibody(2C5)Featured Anti-DDDDK Tag Rabbit Polyclonal Antibody, HRP Conjugated, specially designed for your Western blot analysis.                                                        |
| <b>Storage and Stability:</b> | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. |