

Antibodies Product Datasheet

Product Name:	GAPDH Monoclonal Antibody(2B8), AbFluor™ 405 Conjugated
Catalog Number:	N1105
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Sheep,Insect,Yeast
Applications:	WB,IF/ICC,IHC-p
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC: 1:200, IF 1:200.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG1
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	GAPDH
Protein Name:	Glyceraldehyde-3-phosphate dehydrogenase
Other Name:	GAPDH
Human Swiss Prot No.:	P04406
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	GAPDH Monoclonal Antibody(2B8) AbFluor™ 405 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.