

Antibodies Product Datasheet

Product Name:	GAPDH Monoclonal Antibody(2B8), Biotin Conjugated
Catalog Number:	N1111
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Sheep,Insect,Yeast
Applications:	WB,IF/ICC,IHC-p
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB: 1:5000, IHC: 1:200.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG1
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	GAPDH
Protein Name:	Glyceraldehyde-3-phosphate dehydrogenase
Other Name:	GAPDH
Human Swiss Prot No.:	P04406
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	GAPDH Monoclonal Antibody(2B8) Biotin conjugated specially designed for your WB or IHC analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.