

Antibodies Product Datasheet

Product Name:	His-Tag Monoclonal Antibody(4E6), HRP Conjugated
Catalog Number:	N1154
Reactivity:	Species independent
Applications:	WB,IP,IF/ICC
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:3000-10000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	His-Tag
Protein Name:	His Tag
Other Name:	His
Human Swiss Prot No.:	
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	His-Tag Monoclonal Antibody(4E6) specially designed for your tagged protein quantification and localization.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.