

Antibodies Product Datasheet

Product Name:	Myc-Tag Monoclonal Antibody(3E8), HRP Conjugated
Catalog Number:	N1155
Reactivity:	Species independent
Applications:	WB,IP,IF/ICC
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:1000-3000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	Myc Tag
Protein Name:	Myc-Tag
Other Name:	Myc
Human Swiss Prot No.:	
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	Myc-Tag Monoclonal Antibody(3E8) specially designed for your tagged protein quantification and localization.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.