

Antibodies Product Datasheet

Product Name:	Myc-Tag Monoclonal Antibody(3E8), AbFluor™ 555 Conjugated
Catalog Number:	N1177
Reactivity:	Species independent
Applications:	WB,IP,IF/ICC
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:1000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	Myc Tag
Protein Name:	Myc-Tag
Other Name:	Myc
Human Swiss Prot No.:	
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	Myc-Tag Monoclonal Antibody(3E8) AbFluor™ 555 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.