

Antibodies Product Datasheet

Product Name:	TBP/TATA Binding Protein Monoclonal Antibody(4H2),HRP Conjugated
Catalog Number:	N1206
Reactivity:	Human,Rat,Mouse
Applications:	WB
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:1000-2000.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	TBP
Protein Name:	TATA-box-binding protein (TATA sequence-binding protein) (TATA-binding factor) (TATA-box factor) (Transcription initiation factor TFIID TBP subunit)
Other Name:	TBP
Human Swiss Prot No.:	P20226
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	TBP/TATA Binding Protein Monoclonal Antibody(4H2) HRP Conjugated, specially designed for your Western blot analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.