

Antibodies Product Datasheet

Product Name:	β -Tubulin Monoclonal Antibody(5G3), AbFluor™ 350 Conjugated
Catalog Number:	N1253
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Sheep,Insect,Yeast
Applications:	WB,IF/ICC,IHC-p
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC: 1:200, IF 1:200.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG1
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	TUBB3
Protein Name:	Tubulin beta-3 chain
Other Name:	TUBB3
Human Swiss Prot No.:	Q13509
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	β -Tubulin Monoclonal Antibody(5G3) AbFluor™ 350 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.