

Antibodies Product Datasheet

Product Name:	His-Tag Monoclonal Antibody(4E6)
Catalog Number:	N1268
Reactivity:	Species independent
Applications:	WB,IP,IF/ICC
Dilutions:	WB: 1:3000 IP: 1:200 IF: 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic peptide:HHHHHH conjugated to KLH.
Concentration:	
Gene Name:	His-Tag
Protein Name:	His Tag
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects C-terminal, internal, and N-terminal His-tag fusion proteins.
Storage and Stability:	-20°C/1 year