

Antibodies Product Datasheet

Product Name:	GST-Tag Monoclonal Antibody(1B10)
Catalog Number:	N1271
Reactivity:	Species independent
Applications:	WB,IP
Dilutions:	WB: 1:5000,IP 1:200-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of GST-Tag
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects C-terminal, internal, and N-terminal GST fusion proteins.
Storage and Stability:	-20°C/1 year