

Antibodies Product Datasheet

Product Name:	Lamin B1 Monoclonal Antibody(7C11)
Catalog Number:	N1300
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	WB: 1:2000-5000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Lamin-B1
Concentration:	
Gene Name:	LMNB1
Protein Name:	Lamin-B1
Other Name:	LMNB1; LMN2; LMNB; Lamin-B1
Human Swiss Prot No.:	P20700
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Lamin B1 protein.
Storage and Stability:	-20°C/1 year