

Antibodies Product Datasheet

Product Name:	β I tubulin Monoclonal Antibody(3F7)
Catalog Number:	N1301
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:500-10000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of β I tubulin
Concentration:	
Gene Name:	TUBB1
Protein Name:	Tubulin beta-1 chain
Other Name:	TUBB1; Tubulin beta-1 chain; TUBB2A; TUBB2; Tubulin beta-2A chain; TUBB2B; Tubulin beta-2B chain; TUBB4B; TUBB2C; Tubulin beta-4B chain; Tubulin beta-2 chain; Tubulin beta-2C chain; TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-II
Human Swiss Prot No.:	Q9H4B7
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenousβ I tubulin proteins.
Storage and Stability:	-20°C/1 year