

Antibodies Product Datasheet

Product Name:	HSP70 Monoclonal Antibody(3G10)
Catalog Number:	N1306
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:1000-2000 IF: 1:100-200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of HSP70
Concentration:	
Gene Name:	HSPA1L/HSPA1A
Protein Name:	Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B
Other Name:	
Human Swiss Prot No.:	P34931/P08107
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous HSP70 proteins.
Storage and Stability:	-20°C/1 year