

Antibodies Product Datasheet

Product Name:	HER2 Monoclonal Antibody(11H9)
Catalog Number:	N1309
Reactivity:	Human,Mouse,Rat
Applications:	WB,IF/ICC,IHC-p
Dilutions:	WB: 1:2000-4000 IHC: 1:200 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of HER2
Concentration:	
Gene Name:	ERBB2
Protein Name:	Receptor tyrosine-protein kinase erbB-2
Other Name:	ERBB2; HER2; MLN19; NEU; NGL; Receptor tyrosine-protein kinase erbB-2; Metastatic lymph node gene 19 protein; MLN 19; Proto-oncogene Neu; Proto-oncogene c-ErbB-2; Tyrosine kinase-type cell surface receptor HER2; p185erbB2; CD340
Human Swiss Prot No.:	P04626
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous ErbB-2/HER-2 proteins.
Storage and Stability:	-20°C/1 year