

Antibodies Product Datasheet

Product Name:	CD45 Monoclonal Antibody(12A9)
Catalog Number:	N1310
Reactivity:	Human
Applications:	IF/ICC, WB, IHC-p
Dilutions:	IF: 1:50-200 WB: 1:2000 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CD45
Concentration:	
Gene Name:	PTPRC
Protein Name:	Receptor-type tyrosine-protein phosphatase C
Other Name:	PTPRC; CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD45
Human Swiss Prot No.:	P08575
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CD45 proteins.
Storage and Stability:	-20°C/1 year